



SLEEP DISORDERS CENTER OF VIRGINIA

Say Goodnight Virginia

Policies & Procedures

Equipment Maintenance

ALL Locations

Initial Procedure Date: 07/1998

Review Date: 03/2010

Approved by: Kathe G. Henke, Ph.D.

This policy is to ensure the proper and safe running of all medical equipment used for sleep monitoring.

A contract has been established with “Biomedical Equipment & Engineering Service Co” for annual electrical safety checks of the following equipment:

Sandman collection system Chassis leakage current (amperage should be less than 100 microamps)

1. All pulse oximeters
2. All Cpap/Bipap units
3. All portable Cpap (Auto-sets) units

Collection PC's are checked by Joe Walton (Chassis leakage current less than 300 microamps)

Technicians are also responsible for the following

1. Monthly visual inspection of equipment

Note: All equipment maintenance reports for all equipment checked will be kept at the Glenside Lab. Anna Rodriguez is the contact person for all equipment maintenance.